

The EZ-Glide Stair Chair with POWERTraxx provides maximum safety and power when transporting a patient up or down stairs. Its electronically-controlled motor, direct drive transmission, and powered tracks push and pull patient loads of up to 500 pounds. The chair and track system carry the weight, increasing patient safety and reducing operator fatigue and risk of injury.



- Maximum safety when going up and down stairs
- Easy to operate with paddlestyle wireless controller
- Built-in electronic speed control limits track speed during descent to a safe, manageable level
- Electronic brake automatically locks tracks during stops on ascent or descent
- 5-position lift bar and handles reduce the strain of lifting and provide maximum control for operators of various heights (wireless motor controls work in any handle position)
- SAE J3043 compliant when combined with the Secure Lock+



Motorized tracks power chair up and down stairs



Extending footrest provides added support and comfort for patients



# **EZ**-Glide<sup>®</sup> with **POWER**Traxx<sup>®</sup>

### Specifications



HEIGHT (MIN)	37.5 in   953 mm
HEIGHT (MAX)	63.5 in   1613 mm
WIDTH (OVERALL)	20.313 in   516 mm
WIDTH (SEAT)	16.5 in   419 mm
DEPTH (FRONT TO BACK)	51 in   1295 mm
LENGTH (FOLDED)	37.5 in   953 mm
WIDTH (FOLDED)	20.313 in   516 mm
DEPTH (FOLDED)	10 in   250 mm
WEIGHT	56.6 lb   25.67 kg
LOAD LIMIT	500 lb   227 kg
POWER SUPPLY	28V Lithium Ion Battery

## Features

- Rechargeable Milwaukee® lithium-ion battery (included) is easy to remove, commercially available, and runs up to 20 flights of stairs per charge
- All-metal frame construction provides durability and long life
- Smooth 4-inch front swivel wheels and extra tall
  6-inch rear wheels do not collect debris and provide maneuverability in confined spaces
- One inch track-to-ground clearance provides smooth rolling over carpet and rough surfaces, as well as easier descent down stairs

- Extending footrest and lower track angle increase patient sense of security
- Patient seating surface is large and has removable ABS plastic seat panels
- Positive locking hinge mechanism ensures chair is locked in place whether deployed for use or folded for carrying and storage
- The POWERTraxx Conversion Kit can retrofit powered tracks to previous models of the EZ-Glide Stair Chair (No. 0832195)

If you are seeking an **evacuation chair**, rather than a stair chair designed for use by EMS professionals, ask your FERNO sales representative about FERNO's line of evacuation chairs and solutions.

# Accessories

#### PATIENT HANDLING ACCESSORIES



**EZ-Glide Headrest** No. 0822430



EZ-Glide IV Pole No. 0821976



Oxy-Clip™ O² Bottle Holder No. 0821977



Kwik Klip Restraint System No. 0314002



Battery Charger No. 0832199 Battery

EZ-Glide Rear Handles

**EZ-Glide Headrest Pouch** 

No. 0821975

(pouch only)

No. 0314117



**STORAGE** 

**Secure Lock™+** (SAE J3403 Compliant) No. 600251006



**EZ-Glide Secure Mount** No. 0822072

# **⊖FERNO**°

FERNOEMS.COM | 877.733.0911